



BioControl Assurance Gold - AOAC 999.08

SCOPE

This method is applicable to all foods when modified for raw meat or heavily contaminated samples.

PRINCIPLES

The detection of *Salmonella* spp. in raw meat samples is broken down into stages as follows:

- **Pre-enrichment in non-selective liquid medium**

A 1:10 dilution of the sample must be pre-enriched in buffered peptone water supplemented with novobiocin¹ at 35-37°C for 18 to 26 h. Buffered peptone water should be warmed to room temperature or to 36 °C for large volumes. For carcass sponges, buffered peptone supplemented with 0.1% novobiocin solution is added to the moistened sponge to bring the total volume to 60-100 ml and the sample incubated at 35-37°C for 18 to 26 h. In the case of sponges BPW need not be warmed to room temperature before being used to re-hydrate the sponge, for all subsequent additions BPW should be warmed to room temperature.
- **Selective enrichment**

Culture from the pre-enrichment broth is inoculated into RV broth (0.1 ml in 10 mL) and TT broth (1 mL in 10 mL). Selective enrichment broths are incubated at 42 ± 0.5°C for 5 to 8 h.
- **Post enrichment**

Cultures from selective liquid media are combined and inoculate into pre-warmed TSB+n² and incubate at 42 ± 0.5°C for 16-20 h.
- **Enzyme immunoassay**

Follow the manufacturer's instructions retaining TSB+n for confirmation of presumptive positive results.
- **Cultural confirmation**

Presumptive positives can be confirmed from the retained TSB+n by streaking onto XLD, HE and BS agar³. Confirmation carried out at an 'off-site' laboratory must be from retained BPW enrichment. Typical colonies are confirmed as *Salmonella* using biochemical and serological tests as outlined in AS 5013.10.

¹ Suspend 0.1g of novobiocin sodium salt in 100 mL of purified water. Filter sterilise (0.2 µm). Solution is stable up to 60-days in a dark bottle at 2-8 °C. Added at the rate of 4mL/225mL of BPW.

² Tryptone soy broth + 0.1% novobiocin (novobiocin added after autoclaving)

³ Xylose Lysine Deoxycholate (XLD), Hektoen enteric (HE), Bismuth Sulphite (BS)

CHECKLIST

Pre-enrichment	Is the buffered peptone water warmed to room temperature (to 36°C for large quantities)?	_____
	Is novobiocin added to BPW?	_____
	Is the correct amount of enrichment broth used for the weight of sample analysed?	_____
	What volume is used for carcass swabs?	_____
	Is a positive control run with each batch of samples analysed?	_____
	Are reference cultures inoculated into primary enrichment broth at a level of 10 to 100 cells?	_____
	Is pre-enrichment done at 36 ± 1°C for 18-26 h?	_____
Selective-enrichment	Are selective enrichment broths incubated at the appropriate temperature?	_____
	Is preparation of TT completed on the day of use?	_____
Post-enrichment	Are selective enrichment cultures combined for analysis?	_____
	Is TSB+n broth incubated at 42 ± 0.5°C for 6-7 h?	_____
	Is TSB+n broth retained for confirmation of presumptive positive samples?	_____
Enzyme immunoassay	Are the manufacturer's instructions available?	_____
	Are reagents stored at 2-8°C?	_____
	Is incubation carried out for 30 min at 35-37°C?	_____
Cultural confirmation	Is <i>Salmonella</i> isolated in-house from TSB+n broth?	_____
	Are XLD, HE and BS agars used for confirmation?	_____
(if applicable)	If an external laboratory is used is it department approved?	_____
	BPW should be supplied to off-site laboratories for confirmation following AS 5013.10	_____
	Are <i>Salmonella</i> confirmed using AS 5013.10 (with regard to biochemical and serological tests)?	_____